



MultiTask Surface Cleaner - MultiClean™ Presaturated Wipes



- Removes solder paste, flux residues, finger prints, uncured epoxy residues, light oils and grease
- Most economical
- Water and alcohol mix
- ESD-safe
- No rip. No lint. Heavy-duty wipe.
- Unique “slam-shut” package for extra long service life

Product Information

The MultiTask Surface Cleaner - *MultiClean*™ presaturated wipes provide fast and safe cleaning in almost any production environment. Using a strong, shred-proof disposable wipe saturated with high purity isopropyl alcohol and deionized water, it is well suited for cleaning medical devices, electronics and any surface that requires hygienic grade cleaning. Versatile, it cleans most solder pastes and fluxes found in electronics production. It also cleans organic and ionic residues, light oils and other contaminants commonly found in the manufacture of medical devices. With excellent materials compatibility it is both safe and easy to use.

Money-saving refills are available from MicroCare (image to the right). These contain the same wipes and fluid as the tubs above, but reduce costs and help protect the environment by encouraging the re-use of packaging.

AVAILABILITY: Stocked locally by MicroCare distributors world-wide. Contact MicroCare for a distributor near you. For more information contact MicroCare directly or visit www.MicroCare.com.

MicroCare Corporation

595 John Downey Drive
New Britain, CT 06051 USA
Tel: 860.827.0626
Fax: 860.827.8105
Toll-Free: 800.638.0125
Email: TechSupport@MicroCare.com

MicroCare Europe Bvba

Havendoklaan 19, Cargovil
Vilvoorde, B-1804 Belgium
Tel: +32 2 251 95 05
Fax: +32 2 400 96 39
Email: EuroSales@MicroCare.com

Technical Details and Cleaning Data

Boiling Point	82°C / 180°F
Chemical Family	Deionized Water & Alcohol
Cleaning Strength (KB value)	25 (est.)
Evaporation Rate	Fast
Percent Volatile	100%
Odor	Slight, Alcohol
Specific Gravity	0.87

Health, Safety and Environmental Data

Flammable	Yes
Flashpoint	12°C / 53°F
NFPA	Health: 1 Fire: 3 Instability: 1
Organic Content (g/L)	350
RoHS and WEEE Complaint	Yes
REACH	Yes*
US EPA SNAP-Approved	Yes

GHS / HCS Pictogram



Packaging

Presaturated Wipe Container	100 - 5" x 8.5" wipes
443 mls. fluid content	MCC-MLCW
Presaturated Wipe Refills	100 - 5" x 8.5" wipes
443 mls. fluid content	MCC-MLCWR
MiniCube.	6 lbs / 2.72 k. MCC-MLCG
Cubitainer.	32 lb / 14.5 k. MCC-MLCP
Drum.	380 lb / 172.4 k. MCC-MLCD



Discover Perfectly Clean
www.MicroCare.com